APPLICATION EXAMPLES OF THE *KRi* Inverted Pendulum PP-300

Application Note Ref: IP-300 Date: March 15, 1996 Prepared by: KV Ling, NTU & KW Lim, UNSW for KentRidge Instruments Pte Ltd

Application Examples of the *KRi* Inverted Pendulum PP-300

IMPORTANT NOTICES

KentRidge Instruments (S) Pte. Ltd. reserves the right to make changes to the information contained in this manual without notice. KentRidge Instruments Pte. Ltd. assumes no responsibility for any errors or consequential damages that may result from the use or misinterpretation of any information contained in this manual.

Acknowledgements

The Inverted Pendulum **PP-300** is designed and manufactured by KentRidge Instruments Pte. Ltd. in collaboration with the Nanyang Technological University, Singapore.

Copyright 1996-1999 KentRidge Instruments Pte Ltd

All Rights Reserved. No part of this document may be reproduced, transcribed, stored in any retrieval system, translated into any language, or transmitted in any form or by any means, electronic, mechanical, photocopying, recording or otherwise without prior written permission from KentRidge Instruments Pte. Ltd.

Document File: IP300A.DOC Printed in Singapore

TABLE OF CONTENTS

Table of Contents	v
Table of Figures	v
Overview	7
Description of Apparatus	7
Example 1 Simulating a Crane at a Construction Site	8
Example 2 Balancing the Inverted Pendulum	9
Example 3 Swing Up the Inverted Pendulum	9
Example 4 Swinging Up And Balance the Inverted Pendulum at a Specific Orientation	0
References 1	0
Other KRi Products 1	1

TABLE OF FIGURES

Figure 1 Inverted Pendulum Model PP-300	. 7
Figure 2 Using the PP-300 to simulate a crane	. 8
Figure 3 Balancing the pendulum	. 9
Figure 4 Swing up and balancing the pendulum	. 9
Figure 5 Balancing the pendulum at specific orientation	10



APPLICATION EXAMPLES FOR THE KRi INVERTED

PENDULUM PP-300

KV Ling and KW Lim

OVERVIEW

This note has been prepared to guide the user on applications of the *KRi* Inverted Pendulum PP-300. It describes four example applications of the apparatus for teaching or research in control systems.

The four examples of applications are (in order of increasing difficulty):

- 1. Simulating a crane at a construction site,
- 2. Balancing the inverted pendulum,
- 3. Swing up and balancing the inverted pendulum, and
- 4. Swing up and balancing the inverted pendulum at a specific orientation.

DESCRIPTION OF APPARATUS

The Inverted KRi Pendulum PP-300 consists of: (1) a short arm which rotates in a horizontal plane, driven by a DC motor mounted vertically; and, (2) a pendulum rod at one end of the arm which can rotate 360° in a vertical plane orthogonal to that of the arm. The angular position and speed of the DC motor are measured by an optical encoder (with 1000 lines resolution) and tachometer. The pendulum itself is attached to the rotating arm via a servo potentiometer which measures the angular position of the pendulum in the vertical plane. A proprietary mechanism is used to provide the electrical connections



Figure 1 Inverted Pendulum Model PP-300

between the servo potentiometer and terminal blocks on the stationary base.

The *KRi* Inverted Pendulum Apparatus PP-300 is a <u>low cost pilot plant designed</u> <u>for teaching advanced control systems theory in the laboratory</u>. Compared to conventional linear pendulums, it occupies very little bench space. The motor and

pendulum also rotates freely (360° clockwise or counter-clockwise) with no mechanical restrictions. The PP-300 incorporates mechanical safety guards.

This apparatus is <u>useful to illustrate how to design control algorithms - ranging</u> <u>from classical controllers to modern fuzzy-logic or neural controllers - for balancing</u> <u>a pendulum</u>. For example, in the swing up mode, open-loop schemes can be compared with feedback-based optimal controllers. In the balancing mode, it allows demonstration of phase-lead compensators, fuzzy-logic controllers, neural controllers, or variable-structure controllers. It also demonstrates the use of linear models for a non-linear system and requires the realization of mode switching in controllers from swing up to balance. Thus one simple and compact apparatus successfully provides a point of comparison amongst the rich variety of classical and modern algorithms.

EXAMPLE 1 SIMULATING A CRANE AT A CONSTRUCTION SITE

With the pendulum at rest in its natural stable position, the apparatus can be used to emulate the behaviour of a crane at a construction site. Using the voltage to the motor amplifier as the manipulated variable, the arm can be swung from one position to another. This would resemble the crane operator moving a load (in this case, the pendulum and the mass at its tip).



The control objective here is to move the load in the shortest possible time while keeping the load as stable as possible (minimum perturbations on angle α of the load), see Figure 2 on the left.

With this setup, it is possible to study control techniques ranging from simple single loop PID control (with the angular position θ of the arm as the measured variable) and inner loop feedback (with the angular position α of the pendulum as an additional measured variable) to state feedback/state observer or fuzzy/neutral techniques.

An important pedagogical feature of this application is that whilst a linearised model can be developed when operated cautiously with "small perturbations", the plant quickly becomes nonlinear for larger angular displacements. This provides a simple but realistic and repeatable environment for the evaluation of controllers.

Figure 2 Using the PP-300 apparatus to simulate a crane

EXAMPLE 2 BALANCING THE INVERTED PENDULUM

In this application, the pendulum is first positioned upright manually, that is, in a position of unstable equilibrium. The controller is then switched in to balance the pendulum and to maintain this balance in the presence of disturbances, see Figure 3. A simple disturbance may be a light tap on the balanced pendulum. A complex disturbance may be gusts of wind (use a fan!).

This setup can be used to study the control of open loop unstable system. It is a demonstration of the stabilizing benefits of feedback control. A range of control techniques ranging from the simple phase advance compensator to neural net controllers can be applied. Some of these are described in the references listed at the end of this note.



EXAMPLE 3 SWING UP THE INVERTED PENDULUM

Figure 3 Balancing the pendulum

In this application, the pendulum is

initially at rest in a stable, hanging down position. The control objective is to swing the pendulum upright and to subsequently maintain it at this balanced position.



Figure 4 Swing up and balancing the pendulum

The control realization requires mode switching in the controller from swing up to

balance. A variety of performance criteria can be applied to the swing up problem. These may be to achieve swing up in minimum time using minimum energy. Only the rotary arm position angle may be used or both angle and velocity measurements may be utilized This setup allows the demonstration of both open loop optimal control or nonlinear feedback control such as variable structure control.

EXAMPLE 4 SWINGING UP AND BALANCE THE INVERTED PENDULUM AT A SPECIFIC ORIENTATION

In this application, Example 3 is extended so that the pendulum is required to swing up and balance at a specific orientation, measured by the angular position θ_b of the horizontal rotary arm.

The problem posed makes for an interesting comparison of strategies. A minimum time performance criteria for example may compare



Figure 5 Balancing pendulum at specific orientation

a strategy which first swings the pendulum upright and then moves the rotary arm to the desired position θ_b (if necessary) while keeping the pendulum balanced with one which moves the arm to the desired position before swing up and balance. Other strategies may allow both performance objectives to be attempted simultaneously. This setup will be an interesting test-bed for the comparison of controllers with multi-objective performance criteria but only one manipulated variable.

REFERENCES

- 1. Furuta K., M. Yamakita M and S Kobayashi, "*Swing-up control of inverted pendulum using pseudo-state feedback*", Proc Instn Mech Engrs Vol.206, No.14, Pt I Journal of Systems and Control Engineering, p.263-269.
- 2. Furuta K., M. Yamakita, S. Kobayashi and M. Nishimura, "A New Inverted Pendulum apparatus for education in Advances in Control Education" IFAC 1991 pp 191-196
- 3. Nelson J. and L. Gordon Kraft. "*Real-time control of an inverted pendulum system using complementary neural network and optimal techniques*", Proc American Control Conference, Baltimore, Maryland, June 1994, p.2553-2554.
- 4. Wiklund M., A. Kristenson, K.J. Astrom, "A new strategy for swinging up and control of an inverted pendulum", IFAC World Congress 1993, Vol. IX.
- Yamakita M, and K. Furuta, "VSS Adaptive Control based on nonlinear model for TITech Pendulum", Proceedings of the 1992 IEEE International Conference on Industrial Electronics, (IECON 92), San Diego, USA ISBN 0-7803-0582-5/1992 IEEE, pp 1488-1493

OTHER KRi PRODUCTS

KentRidge Instruments Pte Ltd offers a family of control apparatus or equipment for teaching and research in control engineering:

- Coupled-Tank Control Apparatus PP-100
- Fan & Plate Control Apparatus PP-200
- FlexiDrive PP-400
- Dual Process Simulator KI-101
- Mixed Signal Test Unit TU-100
- Controller Boards UC96

For more information, please contact:

